

SYBR Green I for PCR

Technical literature is available at: www.mesgenbio.com. E-mail MesGen Technical Services if you have questions on use of this system: tech@mesgenbio.com



Catalog Number : MPY6155

Lot Number : Refer to vial

Packaging Size : 1000µL □ 5000µL □

Ex (nm) : 497 Em (nm) : 525 Solvent : DMSO

Background

SYBR Green I, is a very sensitive dye for the detection of double stranded DNA (dsDNA), So it has been widely used in non-specific detection of amplification in realtime qPCR experiments. The double-strand DNA-specific SYBR Green I fluorescent reporter offers distinct advantages. SYBR Green I dye is inexpensive, easy to use, and sensitive. Well-designed primers must be used in SYBR Green quantitative RT-PCR reactions because SYBR Green I dye will detect nonspecific products, resulting in an overestimation of the target concentration.

Storage instruction

2-8°C & Protect from light

1. The following table 1 is our SYBR Green I RT-PCR Reagents Kit recipe for reference only. Please optimize it by yourself.

Reagent	Final concentration in the mix
dNTP	0.25mM
Tween20	1%
BSA	0.1% vol
Tris(pH8.4)	50mM
Hot-start Taq DNA polymerase	1.25 u per reaction
NH ₄ Cl	10mM
KCl	20mM
MaCl ₂	2.5mM
Sybr Green	2X

2. On ice, prepare a 2x master mix containing no DNA, by mixing the components in the following order : water, DMSO, Taq polymerase buffer, dNTPs, MgCl₂, Sybr Green, Taq polymerase.

3. Transfer 2x SYBR Green PCR Master Mix to tubes or plates. Add RNase-free water. Mix the individual solutions.

4. Prepare a reaction mix according to Table 2.

Due to the hot start, it is not necessary to keep samples on ice during reaction setup or while programming the real-time cycler.

PCR Reaction Setup: **Table 2**

DNA	Template DNA (<500 ng/reaction)
SYBR Mix	25.0 µL
Primer 1	2ul (5uM)
Primer 2	2ul (5uM)
dd H ₂ O	Added to 50.0 µL

Program your real-time cycler according to the program outlined in **Table 3**

PCR initial activation step	Heating cycling Number of cycles 40
95 °C 5min	96 °C 10s 60 °C 30s

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.