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Alcian Blue Staining Kit

Do not eat Store at +2 to +8° C & in the dark.

Cat.No. MAS0981

Size :2×100ml

Technical literature is available at: www.mesgenbio.com.E-mail MesGen Technical Services if you have questions on use of this system: tech@mesgenbio.com

Description

Alcian Blue is a dye used to detect cell chondrogenesis and it stains the sulfated proteoglycan in cartilage tissue. The Alcian Blue Staining Kit contains 1% Alcian Blue Solution, 0.1% Nuclear Fast Red in convenient, ready-to-use solutions. At pH 2.5, Alcian Blue stains acid mucopolysaccharides and shows a blue color, while Nuclear Fast Red stains nuclei pink to red and the cytoplasm pale pink.

Composition

| | |
|-------------------------|-------|
| 1% Alcian Blue Solution | 100ml |
| 0.1% Nuclear Fast Red | 100ml |

Procedure

1. Fix the cells with 10% buffered formalin for 10 min at room temperature.
2. Rinse the slides in 1X PBS, stain with Alcian Blue solution for 30 min at room temperature.
3. Wash the slides in running tap water for 5 min, rinse in pure water.
4. Counterstain in Nuclear Fast Red solution for 5 min at room temperature.
5. Wash slides three times in pure water.
6. Gradient elution using ethanol, then Xylene transparent 5-10min, and sealed by Neutral gum.

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
For Research Use Only. Not For Use In Diagnostic Procedures.