For Research Use Only. Not For Use In Diagnostic Procedures

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Cat.No.MG1925

Size : 10 Mini-Gels

Technical literature is available at: www.mesgenbio.com.E-mail MesGen Technical Services if you have questions on use of this system: tech@mesgenbio.com

Description

This kit method depends on the hydrolysis of the phosphoprotein phosphoester linkage using sodium hydroxide in the presence of calcium ions. The gel containing the newly formed insoluble calcium phosphate is then treated with ammonium molybdate in dilute nitric acid. The resultant insoluble nitrophospho-molybdate complex is stained with Methyl Green. After destaining, the phosphoproteins are colored green to green-blue. The detection limit is in the nanogram range, but depends on the degree of phosphorylation of the protein. This method will detect the phosphoproteins phosphotyrosine and β -casein in the 40-80 ng/band and 80-160 ng/band range, respectively. The method presented here is for staining minigels. Volumes will need to be increased for larger gels.

Procedure

1. Place the gel in a dish containing 50 ml of deionized H₂O. Agitate the gel on a reciprocal or orbital shaker table for 10 minutes.
2. Move the gel to a dish containing 25 ml of **Solution 1**, and agitate the gel on a reciprocal or orbital shaker table for 15 minutes.
3. Place the gel in 25 ml of **Solution 2**, and agitate the gel on a reciprocal or orbital shaker table for 30 minutes.
4. Rinse the gel rapidly with deionized H₂O to remove the calcium chloride adhering to the surface of the gel.

5. Transfer the gel to 25 ml of **Solution 3**, cover the tray with a lid, and incubate it for 20 minutes at 65°C.
6. Place the gel in 25 ml of **Solution 4**, and agitate the gel on a reciprocal or orbital shaker table for 10 minutes.
7. Repeat Step 6.
8. Transfer the gel to 25 ml of **Solution 5**, and agitate the gel on a reciprocal or orbital shaker table for 20 minutes.
9. Place the gel in 25 ml of **Solution 6**, and agitate the gel on a reciprocal or orbital shaker table for 20 minutes.
10. To destain, transfer the gel to 25 ml of **Solution 1**, and agitate the gel on a reciprocal or orbital shaker table for 15 minutes.
11. Repeat Step 10. At this time, green spots of phosphoproteins should be seen.
12. To completely destain the gel, place it in 25 ml of **Solution 7**, and agitate the gel on a reciprocal or orbital shaker table overnight. During the destaining, change the acetic acid once, after the solution has become green in color.
13. Following visualization of the phosphoproteins, stain the total proteins on the gel using a Coomassie-Blue-staining method.

Storage condition

2-8°C

For Research Use Only. Not For Use In Diagnostic Procedures.

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组分	体积
Solution 1 6X	125ml
Solution 2 1X	250ml
Solution 3 4X	65ml
Solution 4 10X	25ml
Solution 5 2X	125ml
Solution 6 5X	50ml
Solution 7 10X	25ml

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