

# Goat Anti- Rabbit IgG (H+L)-FITC Conjugated

Technical literature is available at: [www.mesgenbio.com](http://www.mesgenbio.com). E-mail MesGen Technical Services if you have questions on use of this system: [tech@mesgenbio.com](mailto:tech@mesgenbio.com)

**MESGEN**  
INNOVATION BIOTECHNOLOGY

**Catalog Number** : MAN2524

**Lot Number** : Refer to vial

**Packaging Size** : 100 $\mu$ L

## Background

This product has been optimized for use as a secondary antibody in fluorescent applications.

## Conjugation

FITC. Ex: 493nm, Em: 528nm

## Overview

**Description** : Goat polyclonal Secondary Antibody to Rabbit IgG - H&L (FITC)

**Reactivity** : Rabbit

**Formulation** : Liquid

**Tested applications** : FC (Flow Cytometry), IF(Immunofluorescent)

FC: The amount of the reagent is suggested to be used  $\leq 1 \mu\text{g}/10^6$  cells . IF: 1/50 ~1/150

Since applications vary, the appropriate dilutions must be used at an assay dependent concentration.

## Properties

**Immunogen** : Rabbit IgG whole molecule

**Clonality** : Goat Polyclonal

**Isotype** : Goat Polyclonal IgG

## Storage instruction

Store at +4°C. and protected from prolonged exposure to light. Do not freeze.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

## Distributor in China ( Mainland )

**Company** : Shanghai hongsheng Biotechnology Co.,Ltd

**Address** : Room 2-205, Xingzhan building, No.185, Liuying Road, Zhabei District, Shanghai

**Tel** : 86-21-56620378 **Fax** : 86-21-56620378

**Email** : [sales@hongshengbio.com](mailto:sales@hongshengbio.com)

**Website** : [www.hongshengbio.com](http://www.hongshengbio.com)