

# Goat Anti- Mouse IgG (H+L)-HRP Conjugated

Technical literature is available at: [www.mesgenbio.com](http://www.mesgenbio.com). E-mail MesGen Technical Services if you have questions on use of this system: [tech@mesgenbio.com](mailto:tech@mesgenbio.com)

**MESGEN**  
INNOVATION BIOTECHNOLOGY

**Catalog Number** : MAN4002

**Lot Number** : Refer to vial

**Packaging Size** : 100µL □ 1000µL □

## Background

This product has been optimized for use as a secondary antibody in chemiluminescent applications.

## Overview

**Description** : Affinity purified Goat Anti-Mouse IgG antibody is conjugated to horseradish peroxidase.

**Reactivity** : Mouse

**Formulation** : Liquid

**Tested applications** : Western blotting 1:2000-1:8000, Dot blot 1:3000-1:8000, Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## Properties

**Immunogen** : The antiserum was developed in goat using the Mouse IgG as the immunogen.

**Clonality** : Goat Polyclonal

## Storage instruction

Store at -20°C. Stable for one year from the date of shipment.

**Note:** The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

## Distributor in China ( Mainland )

**Company** : Shanghai hongsheng Biotechnology Co.,Ltd

**Address** : Room 2-205, Xingzhan building, No.185, Liuying Road, Zhabei District, Shanghai

**Tel** : 86-21-56620378 **Fax** : 86-21-56620378

**Email** : [sales@hongshengbio.com](mailto:sales@hongshengbio.com)

**Website** : [www.hongshengbio.com](http://www.hongshengbio.com)