

Unstained Protein Ladder (Middle Range)

MESGEN
INNOVATION BIOTECHNOLOGY


Cat.No.PM10120

Lot.No.Refer to Vial

Size : 1×500µl , 2×500µl , 10×500µl

Description

Protein Standard (for SDS-PAGE) is a protein standard specially designed by MesGen Biotech. It consists of seven bands with molecular weight of 10 kD, 20 kD, 30 kD, 40 kD, 60 kD, 80 kD and 120 kD. All bands are highly purified, thus exhibit excellent performance in SDS-PAGE gels.



Separated on 4-20% SDS-PAGE and Coomassie blue R-250 stained

Recommended loading volume

Load 5µl per well for a mini-gel.

Storage condition

This product is stable at -20 °C for 12 months.



For Laboratory Use