

Indocyanine Green

MG87552

Technical literature is available at: www.mesgenbio.com. E-mail MesGen Technical Services if you have questions on use of this system: tech@mesgenbio.com

MESGEN
INNOVATION BIOTECHNOLOGY

Packaging Size : 200mg 1g

Lot Number : Refer to vial

Synonym : Cardiogreen, Foxgreen, IC Green, Indocyanine green

CAS : 3599-32-4


Formula : C₄₃H₄₇N₂NaO₆S₂

Molecular Weight : 774.96

Beilstein Registry Number 4115884

MDL number MFCD00013078

Soluble in water (1 mg/ml), DMSO, methanol, and ethanol (less).



Biochem/physiol Actions

Cardiogreen is a cyanine dye which binds tightly to plasma proteins.

Storage condition

Room temperature

For Research Use Only. Not For Use In Diagnostic Procedures.