

Ingredients	mg/L
<b>INORGANIC SALTS</b>	
Calcium chloride dihydrate	1324.62
Magnesium chloride anhydrous	1068.2
Magnesium sulphate anhydrous	1357.66
Potassium chloride	2240
Sodium dihydrogen phosphate anhydrous	876.92
<b>AMINO ACIDS</b>	
DL-Serine	1100
Glycine	650
L-Alanine	225
L-Arginine hydrochloride	700
L-Asparagine	350
L-Aspartic acid	350
L-Cystine dihydrochloride	25
L-Glutamic acid	600
L-Glutamine	600
L-Histidine hydrochloride	3377.64
L-Isoleucine	50
L-Leucine	75
L-Lysine hydrochloride	625
L-Methionine	50
L-Phenylalanine	150
L-Proline	350
L-Threonine	175
L-Tryptophan	100
L-Tyrosine disodium salt	72
L-Valine	100
β-Alanine	200
<b>VITAMINS</b>	
Choline chloride	0.2
D-Biotin	0.01
D-Calcium-Pantothenate	0.02
Folic acid	0.02
Niacin	0.02
Pyridoxine hydrochloride	0.02
Riboflavin	0.02
Thiamine hydrochloride	0.02
myo-Inositol	0.02
p-Amino benzoic acid (PABA)	0.02
<b>OTHERS</b>	

Alpha-Ketoglutaric acid	370
D(+) Glucose	700
D-Fructose	400
Fumaric acid	55
L-Malic acid	670
Lactalbumin hydrolysate	3330
Succinic acid	60
Sucrose	26680
Yeast extract	3330
pH 6.2	