

Inorganic Salts	
成分	g/L
Calcium Chloride	0.1854
Magnesium Chloride • 6H ₂ O	0.1649
Potassium Chloride	0.4
Potassium Phosphate Dibasic	0.06
Sodium Bicarbonate	1.35
Sodium Chloride	8
Sodium Phosphate Dibasic	0.0479
Amino Acids	
Alanine	0.0356
Arginine • HCl	0.3175
Asparagine • H ₂ O	0.06
L-Glutamine	0.584
L-Aspartic Acid	0.0532
Cystine • 2HCl	0.0782
Glutamic Acid	0.0588
Glycine	0.03
Histidine • HCl	0.1047
Isoleucine	0.13
Leucine	0.13
L-Lysine • HCl	0.1812
Methionine	0.0375
Phenylalanine	0.08
Proline	0.046
Serine	0.042
Threonine	0.12
Tryptophan	0.025
Tyrosine • 2Na	0.1189
Valine	0.115
Vitamins	
D-Pantothenic Acid • ½Ca	0.002
Choline Chloride	0.002
Folic Acid	0.002
myo-Inositol	0.004
Niacinamide	0.002
Pyridoxal • HCl	0.002
Riboflavin	0.0002
Thiamine • HCl	0.002
Other	
Glucose	1
Phenol Red • Na	0.01
Sodium Pyruvate	0.275