

Ingredients	mg/L
<b>INORGANIC SALTS</b>	
Calcium chloride dihydrate	1319.82
Magnesium chloride anhydrous	1068.2
Magnesium sulphate anhydrous	1357.63
Potassium chloride	2870
Sodium bicarbonate	350
Sodium phosphate monobasic	876.92
<b>AMINO ACIDS</b>	
Glycine	650
L-Alanine	225
L-Arginine hydrochloride	700
L-Asparagine	350
L-Aspartic acid	350
L-Cystine hydrochloride	25
L-Glutamic acid	600
L-Glutamine	600
L-Histidine hydrochloride	3086.34
L-Isoleucine	50
L-Leucine	75
L-Lysine hydrochloride	625
L-Methionine	50
L-Phenylalanine	150
L-Proline	350
L-Serine	550
L-Threonine	175
L-Tryptophan	100
L-Tyrosine disodium salt	62.09
L-Valine	100
<b>VITAMINS</b>	
Choline chloride	0.2
D-Biotin	0.01
D-Pantothenic acid (hemicalcium)	0.02
Folic acid	0.02
Niacin	0.02
Pyridoxine hydrochloride	0.02
Riboflavin	0.02
Thiamine hydrochloride	0.52
myo-Inositol	0.02
p-Amino benzoic acid (PABA)	0.02
<b>OTHERS</b>	
D(+) Glucose	1000
Tryptose Broth	2600