

Ingredients	mg/L
INORGANIC SALTS	
Calcium chloride dihydrate	265
Magnesium sulphate anhydrous	96.3
Potassium chloride	400
Sodium chloride	6800
Sodium dihydrogen phosphate anhydrous	122
NaHCO ₃	2200
AMINO ACIDS	
DL-Valine	46.75
Glycine	11.25
L-Alanine	26.75
L-Arginine hydrochloride	167.5
L-Aspartic acid	20
L-Cystine	12
L-Glutamine	292
L-Histidine hydrochloride monohydrate	10.5
L-Isoleucine	26
L-Leucine	26
L-Lysine hydrochloride	37
L-Methionine	15
L-Phenylalanine	16.5
L-Proline	17.25
L-Serine	21
L-Threonine	48
L-Tryptophan	10
L-Tyrosine disodium salt	26.12
VITAMINS	
Choline bitartrate	2.535
D-Pantothenic acid (hemicalcium)	1.2
Folic acid	2.2
Niacinamide	0.6
Pyridoxal hydrochloride	1
Riboflavin	0.375
Thiamine hydrochloride	1.7
myo-Inositol	2.15
OTHERS	
Glucose	1000